

Stance Control Orthotic Knee Joint Now Available.

This knee joint locks at heel strike and releases at heel off.

Due to the complexity of the knee joint, a training program will be required and will be open to ABC-certified or state-licensed practitioners prior to purchasing. John Michael, MEd, CPO, FISPO, FAAOP is the instructor.

Designed for patients with weak quadriceps and varying degrees of knee instability, this lightweight orthosis will allow patients to regain their mobility and assist them in a more energy-efficient ambulation.

For information or course dates, call 1.866.663.3970, or visit our Web site at www.stancecontrol.com.



Polio Patient



Features:

- Lightweight materials
- High durability
- Automatic lock at heel strike
- Free motion during swing phase
- Can be switched from free knee to 180-degree locked extension or stance control mode.



XXXXXXXX

Horton Technology Inc.

Innovations In Orthotics and Prosthetics

1.866.663.3970
www.stancecontrol.com

Circle 94 on Reader Service Card